ANDERSZ www.andersz.com tel. 502 652 932

VPS Impression Materials Ideal Materials For Every Technique & Speed

Highest quality and value priced. Developed with the help of leading clinicians to make consistent, reliable, results easy! Wash materials are super-hydrophilic and give a high-gloss surface and high-tear strength. Precise setting time and snap-set. Placement in mouth triggers an accelerated set. Formulated from the highest quality polymers and fillers.

Available in six color coded viscosities for easy differentiation: Super Flow, Light Body, Monophase, Heavy Flex, Heavy Stiff Bite, Clear Bite & Putty.

The material is free of taste and is suitable for a variety of impression-taking techniques.

- Putty and Super Flow Light Body formulations can be employed for the standard putty-wash technique; the smooth putty resists sticking to fingers and gloves.
- Heavy Stiff Bite is the perfect bite registration material and provides a high degree of accuracy.
- Heavy Flex can be used for closed bite impressions (using Triple Tray impression trays) in combination with Super Flow Light Body for closed-bite impressions.
- The very fast-setting Monophase material also is formulated for Triple Tray impressions and is particularly suitable for those that do not involve hard-to-reach areas (which generally require a super fl ow viscosity).

3 SET TIMES



First Quarter Fast Set Impression Material

First Half

First Quarter

First Quarte

Star VPS

Star VPS

irst Hah

A super fast impression material great for smaller, Triple Tray type impressions.

Star VPS

Star VPS




First Half Medium Set Impression Material

Ideal for those that want speed, but need more time than is offered by the First Quarter series. New Orthodontic Monophase is ideal for single viscosity full arch impressions. Available in 50ml cartridges or 380ml cartridges for Pentamix.



Star VPS Regular Set Impression Material

offers a variety of viscosities in the standard set times, leaving plenty of working time for multiple unit impressions!





ANDERSZ www.andersz.com tel. 502 652 932

First Quarter	Working Time	Total Time To Set	Durometer	Color	2 Pack4	Pack	Econo 24 Padk	Super Ecoro 96 Pack
Light Body	1/2 Minut e	3/4 Minute s	48		89277	89272	89462	90731
Monophase (Chocomint Flavor)			62		90901	90900	90902	90903
Monophase			62		89386	89387	80008-01	90729
Heavy Body			65		90913	90809	90810	90811

First Half	Working Time	Total Time To Set	Durometer	Color	2 Pack	4 Pack	Econo 24 Padk	Super Ecoro 96 Pack
Light Body			48		90467	90468	90469	90733
Light Body			48		91241	91242	91243	91244
Monophase Ortho		21/4 Minute s	62		91122	91123	91124	91125
Heavy Flex			63		90648	90649	90650	90732
Putty-FS			73		93133	(2 400gm Tubs Base & Catalys)		

Star VPS	Working Time	Total Time To Set	Durometer	Color	2 Pack	4 Pack	Econo 24 Padk	Super Ecoro 96 Pack
Ultra Light	21/2 Minute s	4 /2 Minute s	48		N/A	80060-01	80061-01	N/A
Light Body			48		80021-01	80011-01	80009-01	90727
Heavy Flex			63	٩	80023-01	80013-01	80007-01	90728
Ultra Heavy			65N	2	/A	80062-01	80065-01	N/A
Heavy Stiff Bite	20 Seconds	50 Seconds	90		80024-01	80014-01	80006-01	90730
Clear Bite	25 Seconds 2	Minutes	63		80025-01	80015-01	80037-01	N/A
Putty	90 Seconds	4 ½ Minutes	65		80016	(2 400gm Tubs Base & Catalys)		
AlphaBite	30 Seconds	60 Seconds	82	٩	93352	N/A	93354N	/A
Noginate	90 Seconds	2 ½ Minutes	45		93553	N/A	93554N	/A
								-



ZEST DENTAL SOLUTIONS

TEL: 800.262.2310

ANDERSZ www.andersz.com tel. 502 652 932 USA WWW.ZESTDENT.COM